

Product Name:	NSBDR492SSG/0
Description:	<p>Black Diamond Series High Performance Laboratory and Pharmacy Glass Door Refrigerator 46 Cu. Ft. capacity Two swing Stainless Steel glass doors, self-closing, right hinged, not reversible Made in USA High performance refrigeration system Heavy duty construction Digitally commutated fans Dual location NIST compliant temperature sample probes Dual temperature sample display Premium internal LED lighting with custom color selection External LED flashing visual alarm Data logging (temperature history) USB exportable data download Extreme environment operation Hyper-recovery operation Probe access port (1") 10 shelves Independent digital microprocessor temperature controller Adjustable operating temperature range: 1°C – 10°C Touchscreen display with intuitive interface:</p> <ul style="list-style-type: none"> • Password protected settings and data • Display includes active state verification • °C/°F readout switch • Temperature chart (24 hours) - based on dual probe temperature readings • Visual temperature scale bar - based on dual probe temperature readings • Decimal temperature readout • Battery backup for temperature display and logging • Time zone specific clock display • Tech support info <p>Optional smart ThermalView™ security glass on door Hydrocarbon, natural refrigerant (R290) Premium locking system Internal plenum air distribution Low profile commercial grade casters Standard temperature test validation Single point HI and LO temperature history Temperature calibration (offset) Audible alarm Remote alarm contacts Power failure alarm Sensor error alarm Low battery alarm Audible alarm ring back Door alarm Manual alarm for test validation Easily cleaned smooth interior Exterior dimensions: 52 ¾" W x 35 1/4" D x 80 ½" H Cycle defrost Energy Star Certified 115V, 60 Hz, 5.9 Amps, ½ HP UL/C-UL listed 3 Year parts and labor warranty, plus an additional 2 year compressor parts warranty Shipping weight: 710 lbs.</p>