

Product Name:	LF141WWW/0HC
Description:	<p>General Purpose Auto Defrost Hydrocarbon Laboratory Freezer 14 Cu. Ft. Capacity One solid swing door, right hinged, not reversible Mechanical thermostat Adjustable operating temperature range: -15°-25°C Three adjustable shelves Keyed door lock Magnetic door gasket for positive seal Heavy duty door handle Leveling legs Probe access port (3/8") White textured cabinet and door High density urethane foam cabinet and door insulation Hydrocarbon, natural refrigerant (R600a) One year parts and labor warranty; plus additional four year compressor parts warranty Exterior dimensions: 28" W x 26 1/4" D x 67 3/4" H Auto defrost 115V, 60Hz, 5 Amps, 1/4 HP UL/C-UL listed Shipping weight: 185 lbs.</p>