

Product Name:	NSPF201WWW/0MHC
Description:	<p>Nor-Lake Scientific Manual Defrost Hydrocarbon Laboratory &amp; Medical Freezer  20 Cu. Ft. Capacity  One solid swing door, right hinged  Microprocessor temperature controller  Digital temperature display  Audible and visual high and low temperature alarms  Remote alarm contacts  2 probes; 1 in air and 1 in glycol bottle  Adjustable operating temperature range: -15°C – -25°C  Six fixed shelves  Keyed door lock  Magnetic door gasket for positive seal  Leveling legs  Probe access port (3/8")  Powder coat exterior  High density urethane foam cabinet and door insulation HFC-free refrigerant (R290), foam insulation and packaging Two year parts and labor warranty, plus additional three year compressor parts warranty  Exterior dimensions: 30 3/4 "W x 28 1/2" D x 74 5/8"H  <i>Note: This unit must have 4" clearance on sides and back for adequate ventilation</i>  Manual defrost  115V, 60 Hz, 1.96 Amps, 1/2 HP  UL/C-UL listed  ENERGY STAR CERTIFIED  Shipping weight: 254 lbs.</p>