

Product Name:	FF141WWW/0MHC
Description:	<p>White Diamond Series Natural Refrigerant Flammable Storage Freezer 14 Cu. Ft. Capacity One solid swing door, right hinged Microprocessor temperature controller Digital temperature display with sample simulated probe Audible and visual high and low temperature alarms Remote alarm contacts Adjustable operating temperature range: -15°C – -25°C Six fixed shelves Keyed door lock Magnetic door gasket for positive seal Leveling legs Probe access port (3/8") Powder coat exterior High density urethane foam cabinet and door insulation HFC-free refrigerant (R290), foam insulation and packaging Three year parts and labor warranty Exterior dimensions: 23 5/8 "W x 24 3/8" D x 73"H Note: This unit must have 4" clearance on sides and back for adequate ventilation Manual defrost 115V, 60 Hz, 1.96 Amps, ½ HP UL/C-UL listed Shipping weight: 185 lbs.</p>