

Product Name:	FF201WWW/0MHC
Description:	<p>White Diamond Series Natural Refrigerant Flammable Storage Freezer</p> <p>20 Cu. Ft. Capacity</p> <p>One solid swing door, right hinged</p> <p>Microprocessor temperature controller</p> <p>Digital temperature display with sample simulated probe</p> <p>Audible and visual high and low temperature alarms</p> <p>Remote alarm contacts</p> <p>Adjustable operating temperature range: -15°C – -25°C</p> <p>Six fixed shelves</p> <p>Keyed door lock</p> <p>Magnetic door gasket for positive seal</p> <p>Leveling legs</p> <p>Probe access port (3/8")</p> <p>Powder coat exterior</p> <p>High density urethane foam cabinet and door insulation</p> <p>HFC-free refrigerant (R290), foam insulation and packaging</p> <p>Three year parts and labor warranty</p> <p>Exterior dimensions: 30 ¾ "W x 28 ½" D x 74 5/8"H</p> <p><i>Note: This unit must have 4" clearance on sides and back for adequate ventilation</i></p> <p>Manual defrost</p> <p>115V, 60 Hz, 1.96 Amps, ½ HP</p> <p>UL/C-UL listed</p> <p>Unit meets NFPA guidelines 45 and 70</p> <p>Shipping weight: 240 lbs.</p>