

Product Name:	PRF072WWG/0A
Description:	<p>White Diamond Refrigerator & Freezer Combination with Auto Defrost Freezer 7 Cu. Ft. capacity (5.2 Ref./1.3 Freezer) Two exterior swing doors (glass door refrigerator/solid door freezer), right hinged, not reversible Microprocessor temperature controller, refrigerator and freezer Digital temperature display with sample simulated probe, refrigerator and freezer Audible and visual high and low temperature alarms, refrigerator and freezer Remote alarm contacts, refrigerator and freezer Adjustable operating temperature range: Refrigerator: 1°C – 10°C Adjustable operating temperature range: Freezer:-15°C – -25°C Auto defrost cycle every 6 hours, max air temp in empty freezer: +4.34°C; max air temp in loaded freezer: -2.89°C. Max sample temp in empty freezer: -4.76°C; max sample temp in loaded freezer: -9.58°C. 3 adjustable shelves refrigerator/2 adjustable shelves freezer, shelves have guard rails to prevent product from being pushed against back wall Footed bottom shelf refrigerator, and freezer Back wall evaporator cover Forced draft circulation Keyed door locks Magnetic door gasket for positive seal Integrated door handles Leveling legs Stacking kit Two probe access ports (3/8"); one, refrigerator/one, freezer White powder coated exterior HFC-free refrigerant (R600a), foam insulation and packaging Three year parts and labor warranty, plus an additional two year compressor parts warranty Exterior dimensions: 23 3/4"W x 24 1/4"D x 53 1/8"H Auto defrost freezer / cycle defrost refrigerator 115V, 60 Hz, 1.3 Amps refrigerator, 1 Amps freezer, 1/6 HP refrigerator, 1/8 HP freezer UL/C-UL listed Shipping weight: 207 lbs. In stock, with next day shipping *Please note this unit requires 2 outlets. DISCLAIMER: Auto defrost freezers incorporate an electric heater installed on the evaporator that assist in the removal of the frost and</p>

	ice. Although the heater helps speed up the frost and ice removal, it can raise the temperatures within the freezer.
--	--