

Product Name:	PR031WWG/0
Description:	<p>White Diamond Series Undercounter Refrigerator Freestanding  2.5 Cu. Ft. Capacity  One glass swing door, right hinged, not reversible  Microprocessor temperature controller  Adjustable operating temperature range: 1°C – 10°C  Digital temperature display  Audible and visual high and low temperature alarms  2 probes; 1 in air and 1 in sample bottle  Remote alarm contacts  Four adjustable shelves (two large, two small) with guard rail on back to prevent product from being pushed against the back wall  Back wall evaporator cover  Forced draft circulation  Keyed door lock  Magnetic door gasket for positive seal  Leveling legs  Probe access port (3/8")  Powder coat exterior  High density urethane foam cabinet and door insulation  HFC-free refrigerant (R600a), foam insulation and packaging  Three year parts and labor warranty, plus an additional two year compressor parts only warranty  Pyxis®, Omnicell® and AcuDose RX® compatible  Exterior dimensions: 17 ¾ "W x 19 ½ "D x 28 7/8"H  <i>Note: This unit must have 4" clearance on sides and back for adequate ventilation</i>  Cycle defrost  115V, 60 Hz, .90 Amps, 1/8 HP  ADA compliant  UL/C-UL listed  Shipping weight: 97 lbs.</p>