

Product Name:	PR031WWW/0FB
Description:	<p>White Diamond Series Undercounter Built-In Refrigerator  2.5 Cu. Ft. Capacity  One solid swing door, right hinged  Microprocessor temperature controller  Adjustable operating temperature range: 1°C – 10°C  Digital temperature display  Audible and visual high and low temperature alarms  2 probes; 1 in air and 1 in sample bottle  Remote alarm contacts  Four adjustable shelves (two large, two small) with guard rail on back to prevent product from being pushed against the back wall  Back wall evaporator cover  Forced draft circulation  Self-closing door  Keyed door lock  Magnetic door gasket for positive seal  Leveling legs  Probe access port (3/8")  Powder coat exterior  High density urethane foam cabinet and door insulation  HFC-free refrigerant (R600a), foam insulation and packaging  Three year parts and labor warranty, plus an additional two year compressor parts only warranty  Pyxis®, Omnicell® and AcuDose RX® compatible  Exterior dimensions: 17 ¾ "W x 22 3/4 "D x 30 3/4"H  Energy Star Certified for Laboratory Refrigerators and Freezers  Cycle defrost  115V, 60 Hz, .90 Amps, 1/8 HP  Energy Consumption 0.58 kWh/day  ADA compliant  UL/C-UL listed  Shipping weight: 85 lbs.</p>